

Laboratory pH Meter Calibration Log

Quality Control

DATE ____/____/____	CALIBRATED BY	BUFFER LOT NUMBERS	BUFFER EXPIRY DATES
------------------------	---------------	--------------------	---------------------

	METER 1	METER 2	METER 3	METER 4	METER 5
Meter ID / Serial					
Electrode ID					
Buffer 1 (pH 4.01)					
Buffer 2 (pH 7.00)					
Buffer 3 (pH 10.01)					
Slope after calibration					
Offset (mV)					
Temperature					
Electrode condition					
Calibration accepted					
Initials					

Notes

Go digital with Miratag

Skip the paper — fill checklists faster on the Miratag mobile app. Add photos, videos, notes and signatures. Track compliance in real time. Start free at miratag.com

This template is a general-purpose resource, not tailored to any specific jurisdiction. Each organisation must validate compliance with local regulations.



Scan to view